



Savitzky-Golay

Define window size in ms. Full window is $2 \times \text{window size} + 1$

Apply the quadratic ($m=2$) function or quartic ($m=4$) as a sliding window along the time series to be smoothed.
The wider the window the poorer the time resolution



Convolve window with time series.



```

procedure TdlgData3DFilter.FilterDataSavGol(var ProcessedArray: TProcessedPosRecArray); //Savitsky-Golay generalised moving average filter that convolves a polynomial with the signal to smooth (low pass filter)
var n,a: single;
    m4: boolean; // The window size defines the cutoff frequency of the low-pass filter
    procedure AddSample(j,t,dt: integer); // The total length of the filter window must be a positive odd integer less than the length of the TimeSeries
        var w: single; //m2 = quadratic polynomial m4 = quartic polynomial
        begin //all times measured in ticks (100,000 ticks per second)
            w:=(ProcessedArray[j].time - t)/dt; // calculate w (window time) as a fraction of the window size
            //(i.e. to lie in the range 0 to 1)

            if m4 then
                w:=-4.45*sqr(w)+3.83*sqr(sqr(w))+1 //calculate the quartic window weight
            else
                w:=-1.58*sqr(w)+1; // or calculate the quadratic window weight
            a:=a+ProcessedArray[j].a*w; //apply the weight
            n:=n+w;
        end;
    var i,j,t,dt: integer;
    ProcessedArray2: TProcessedPosRecArray;
    begin
        ProcessedArray2:=copy(ProcessedArray);
        dt:=seWidth.Value*TicksPerSec div 1000; //sewidth is window size in ms
        // sewidth of 1ms = 1*100,000/1000 = 100ticks

        m4:=rbM4.Checked;
        for i:=low(ProcessedArray) to high(ProcessedArray) do //ProcessedArray contains signal times series to be filtered
        begin
            n:=1;
            a:=ProcessedArray[i].a; //signal amplitude
            t:=ProcessedArray[i].time; //frame timestamp in ticks (in AAA samples are not forced to be evenly spaced in time although they mostly are).

            j:=i-1;
            while (j >= low(ProcessedArray)) and (ProcessedArray[j].time >= t-dt) do //Filters forward in time. Does not filter samples that are less than the specified window width (in ms) from the sample being smoothed
            begin //i.e. applies the left side of the symmetric window function
                AddSample(j,t,dt);
                dec(j);
            end;

            j:=i+1;
            while (j <= high(ProcessedArray)) and (ProcessedArray[j].time <= t+dt) do //Filters backward in time. Cannot filter samples that are more than the specified window width (in ms) from the sample being smoothed
            begin //i.e. applies right side of the symmetric window function
                AddSample(j,t,dt);
                inc(j);
            end;

            ProcessedArray2[i].a:=a/n; //divide by cumulative weight total so as not to change the amplitude.
        end;

        move(ProcessedArray2[0],ProcessedArray[0],sizeof(ProcessedArray[0])*length(ProcessedArray));

        PosVelAcc(ProcessedArray);
    end;

```